

*Unique, Patented Nasal Powder Spray
Relieves Symptoms **Without Side Effects!***

FastBlockAllergy.com



Exclusively From:

*Terry
Naturally*

STOP ALLERGIES...Fast Relief Begins Here!

FastBlock™ Allergy Relief



WHAT IS FASTBLOCK™ ALLERGY RELIEF?

FastBlock Allergy Relief is unlike any other treatment for allergies. It's not a liquid. It's not an antihistamine or steroid. It's the **WORLD'S ONLY DRUG-FREE POWDER ALLERGY SPRAY**. But unlike liquid nasal sprays or oral medications, FastBlock Allergy Relief doesn't cause drowsiness, irritate the nasal passages, or cause unpleasant side effects.

FastBlock Allergy Relief is clinically proven to deliver fast, effective protection against mild allergies and hay fever caused by airborne allergens such as pollen, dust mites, and animal dander. It's a patented device that distributes a fine cellulose powder into the nasal passages, coating the allergy receptors in the nose where airborne allergies attack.



WHO SHOULD USE IT?

Anyone suffering from hay fever or allergies caused by airborne allergens can use FastBlock Allergy Relief. Whether from seasonal allergies and hay fever or everyday allergies to pets, or dust mites—it will relieve symptoms such as sneezing, sinus congestion, and runny, itchy nose. Because FastBlock Allergy Relief is completely safe, it's suitable for children as young as 18 months along with older adults and pregnant or nursing women. It's easy to use and can be repeated as needed.

WHO WILL IT HELP?

50 MILLION

Nasal Allergies affect over 50 million people in the United States.

Source: Asthma and Allergy Foundation of America

17.6 MILLION

According to the CDC, in 2012 there were 17.6 million people diagnosed with hay fever.

Source: American College of Allergy, Asthma, and Immunology

10 MILLION

Roughly 10 million adults and children who have pets are actually allergic to them.

Source: American College of Allergy, Asthma, and Immunology

13.4 MILLION

More than 13.4 million annual visits to physician offices, hospital outpatient departments and emergency rooms were due to allergic rhinitis.

source: American College of Allergy, Asthma, and Immunology

FREQUENTLY ASKED QUESTIONS...

INDICATIONS

Relieves mild allergic symptoms of hay fever and allergic rhinitis: sneezing, itchy nasal passages, runny nose, and sinus congestion.


HOW DOES IT WORK?

FastBlock Allergy Relief works with the body's own defense mechanisms to strengthen resistance to airborne allergens. The unique delivery system administers a fine powder spray into the nasal passages, allowing the cellulose powder to mix with the natural moisture in your nose to form a protective, gel-like barrier. This barrier prevents airborne allergens from coming into contact with the allergy receptor sites in the nasal passages, thereby avoiding an allergic reaction.

Over 20 positive clinical trials with statistically significant results have been conducted on FastBlock Allergy Relief. It's proven to effectively relieve symptoms such as sneezing, runny, itchy nose, and sinus congestion.



VISIT FASTBLOCKALLERGY.COM
TO VIEW THE VIDEO



FastBlock™ Allergy Relief
patented delivery system

Q. Do I need to worry about interactions with drugs or supplements?

A. No. FastBlock Allergy Relief contains no antihistamines, steroids, or any other form of drug.

Q. When should I use FastBlock Allergy Relief?

A. You can use FastBlock Allergy Relief as a preventive when pollen counts are high, as soon as you notice symptoms, or if you know you'll be in an environment with airborne allergens.

The patented powder and delivery system of FastBlock Allergy Relief has been clinically proven to be effective—without the jitteriness, drowsiness, or other side effects of antihistamines, steroid drugs, and medications.

Q. Why does each spray of FastBlock Allergy Relief smell like peppermint?

A. The light peppermint smell lets you know that you've successfully sprayed it into your nose and that it is going to work to block allergens.

Q. Can I use FastBlock Allergy Relief more than three times a day?

A. Yes. It can be used as often as needed. But chances are you'll need to use it less often as the season goes on because clinical research has shown that FastBlock Allergy Relief becomes more effective with continued use at the same standard dosage throughout the season.

Q. Can children use FastBlock Allergy Relief?

A. Yes. Children 18 months and older can use FastBlock Allergy Relief.




PROVEN SUCCESS OF **FASTBLOCK™** ALLERGY RELIEF IN CLINICAL TRIALS

Significantly Reduces Grass Pollen Allergy Symptoms

RESULTS

- Highly significant reduction in allergy symptoms: sneezing, stuffy or runny nose, irritated eyes
- 87% of participants taking FastBlock Allergy Relief rated their treatment as “good” or “very good” versus only 24% in the placebo group
- Double-blind, placebo-controlled study
- 108 patients (ages 18-40) with a history of seasonal allergies
- 4-week study during grass pollen season
- Participants used three puffs of FastBlock Allergy Relief daily or 3 puffs of a placebo (lactose powder)
- At the end of each day, participants rated the severity of their allergy symptoms
- Results demonstrated a significant reduction in allergy symptoms by the participants using FastBlock Allergy Relief




Åberg, Ospanova, Nikotin, Emberlin, Dahl. *INT Arch Allergy Immunol* 2014;163:313-318.
Sum of nasal symptoms day by day in the respective groups (full analysis set, n = 107).

Reduces Symptoms of House Mite Dust Allergy

RESULTS

- Significant improvement in symptom score for participants using FastBlock Allergy Relief
- Two and a half times more participants saw a significant reduction in their allergy symptoms when using FastBlock Allergy Relief than when using the placebo
- Double-blind, placebo-controlled crossover trial
- 15 participants with known sensitivity to house dust mite allergen
- In phase one of the test, each participant was treated with either FastBlock Allergy Relief or a placebo powder and then exposed to house dust mite allergen; patients were not told which treatment they received
- Allergy symptoms were recorded at regular intervals for 24 hours after exposure to house dust mites
- A week later the test was repeated with the treatment the patient did not receive in phase one
- Use of FastBlock Allergy Relief was associated with a reduction in sneezing, itching and runny nose, and less nasal congestion



Emberlin J, Lewis R. *Curr Med Res Opin* 2007;23(10): 22423-2431.

Significantly Reduces Hay Fever Symptoms in Small Children

RESULTS

- Effectively reduced allergy symptoms in children
 - Particularly effective in reducing sneezing and runny nose
 - No adverse effects
- Double-blind, placebo-controlled study
 - 53 children ages 8–18 years of age
 - 4-week study during springtime allergy season
 - Participants used one puff of FastBlock Allergy Relief three times daily, or placebo
 - A significant reduction in total symptom scores from the nose, especially for runny nose, was experienced by the children using FastBlock Allergy Relief; no significant adverse effects were reported


SYMPTOM	TREATMENT	AVERAGE SYMPTOM RATING
Sneezing	Placebo	2.31
	FastBlock™ Allergy Relief	1.91
Runny Nose	Placebo	2.56
	FastBlock™ Allergy Relief	2.03
Stuffy Nose	Placebo	.42
	FastBlock™ Allergy Relief	2.13
Eye Symptoms	Placebo	2.26
	FastBlock™ Allergy Relief	2.11
Lower Airways	Placebo	.63
	FastBlock™ Allergy Relief	1.47
Sum of All Symptoms	Placebo	11.17
	FastBlock™ Allergy Relief	9.66
Sum of Nasal Symptoms	Placebo	7.29
	FastBlock™ Allergy Relief	6.07

Aberg N, Benson M. *Pediatr Allergy Immunol.* 2011; 22: 594–599

Improves Quality of Life for Hay Fever Sufferers

RESULTS

- Statistically significant reduction in itchy nose, runny nose, sneezing, and nasal congestion
 - Improvement experienced after the first administration of FastBlock Allergy Relief
 - Participants also reported better sleep, better mood, and reduction in allergy symptoms
- 30 people with weed allergy (pollen, ragweed, and plants in the sunflower family)
 - Group 1 received standard allergy treatment
 - Group 2 received standard allergy treatment plus FastBlock Allergy Relief (3 times per day)
 - 4-week study
 - The participants in Group 1 receiving standard allergy treatment had a statistically significant improvement in runny nose and congestion only, and no improvement in quality of life measurements
 - The participants in Group 2, who received standard allergy medication plus FastBlock Allergy Relief, had better improvement of their allergy symptoms as well as improvement in many quality of life measurements, including sleep and mood




Dynamics of allergic rhinitis symptoms among patients who received FastBlock™ Allergy Relief. Penechko Em, Sizyakina LP, *Russian Allergy Journal.* 2011

Blocks Allergens

RESULTS

- Blocked the penetration of dust mite allergens
 - Reduced the body's overall allergen load
- In vitro test of the ability of FastBlock Allergy Relief to block the passage of dust mite allergens
 - Amount of allergen to reach the test solution was less than 20 ng/ml for five hours after exposure
 - Without FastBlock Allergy Relief, the allergen content was almost 100 ng/ml one hour after exposure



Amount of Der p 1 (dust mite allergen) penetrating a 1.5 mm thick HPMC gel layer compared to control (no barrier). Diethart, Emberlin, Lewis. Presented as a Poster at EAACI XXVII Congress, Barcelona, Spain 7-11 June 2008

DETAILS OF CLINICAL TRIALS ON **FASTBLOCK™** ALLERGY RELIEF

STUDY	DESCRIPTION	PARTICIPANTS	RESULTS
Use of Cellulose Powder for the Treatment of Seasonal Allergic Rhinitis <i>Adv Ther.</i> 2003 Jul-Aug;20(4):210-9. Open Clinical Trial	Josling, Steadman. Participants were instructed to use one puff in each nostril of FastBlock Allergy Relief. Drug treatment was allowed if a full hay fever attack occurred. Participants kept a diary documenting the effectiveness of the powder (a score of 5 represents no symptoms and complete control) and average time to relieve symptoms.	102 participants 66 female, 36 male Average age 44	Overall average daily score was reported as 3.85 (on the 5 point scale) for men and women indicating a minimum 77% success rate in chronic hay fever sufferers. Only 12% of volunteers recorded a daily average score under 2.9 revealing a higher score than all pharmaceutical alternatives that were used by participants. Relief was obtained within 0.1-3 hours.
Clinical study of FastBlock Allergy Relief for relief of allergy symptoms including sneezing, runny nose, itchy and watery eyes Presented at the Pan-Hellenic Conference of ENT Specialists on 19th March 2004 in Thessaloniki, Greece. Open Clinical Trial	Vlahtis, K. The product was used once a day, usually in the morning or shortly before the known time of day when symptoms usually appear. One application per nostril. Evaluation at time 0, 3 weeks and 6 weeks. Trial took place from December 2003 till March 2004 to examine participants that were perennial or chronic allergy sufferers rather than seasonal sufferers.	40 participants 24 female, 16 male	All participants were using a pharmaceutical treatment (Decongestant 35%, corticosteroids 42.5%, antihistamines 2.5%, corticosteroid/antihistamine combination 20%) at the beginning of the study. Participants were asked to discontinue use of this medication during the study. After 3 weeks of use, 85% of participants realized improvement in their allergy symptoms. After 6 weeks of use, 90% of participants realized improvement in their symptoms. No reported side effects.
Measure of improvement in nasal mucociliary clearance and PNIFR (peak nasal inspiratory flow rate) in children with allergic rhinitis Presented at World Allergy Congress in Munich, Germany June 2005 Published in <i>Nea Paediatrica Chronica</i> , June 2005, Vol 5 no 2. Open Clinical Trial	Aivazis, Bourli, Maratou et al. Conducted at the University of Thessaloniki, Greece. The mucociliary clearance was determined with use of a non-invasive dye (Edicol Orange 3% + CaHPO ₄ 2H ₂ O 97%) in vivo. The participant's mucociliary clearance was measured prior to using FastBlock Allergy Relief and then 2 days after finishing the 6-week treatment.	100 participants Average age 8.2 Age range 1.5–8.2	There was significant improvement in Nasal Mucociliary Clearance (reduced from 39 minutes to 18.15 minutes) and PNIFR values. The average clearance time was reduced from 55.23 minutes to 21.2 minutes. Out of the 51 children who started the trial with abnormally prolonged clearance, all but 5 children had improved nasal mucous clearance times. No adverse effects were reported.
A double-blind placebo-controlled trial of inert cellulose powder for the relief of symptoms of hay fever in adults <i>Current Medical Research and Opinion</i> 2006;22(2):75-85 Poster presented at World Allergy Congress, Munich 2005	Emberlin, Lewis. A double-blind, placebo-controlled trial in 2004. Participants could take any other medication along with FastBlock Allergy Relief or the placebo.	97 participants 57 female, 40 male Age range 18+	The amount of rescue medication used by the placebo group (overall and in individual categories e.g antihistamines, nasal sprays and eye drops) was significantly greater than that used by the active (FastBlock Allergy Relief) group. No significant difference found in the symptom scores for the two groups. FastBlock Allergy Relief significantly reduced the need for rescue medication. No adverse effects were reported.
Double-blind placebo-controlled cross over trial of inert cellulose powder, by nasal provocation with grass pollen to assess efficacy of the product in controlling symptoms of hay fever Poster presented at EAACI, Vienna June 2006	Emberlin, Lewis. A trial of FastBlock Allergy Relief by nasal provocation with grass pollen (350 grains per cubic meter). At baseline and at regular intervals after challenge scores were taken for 6 symptom categories, measures were taken of nasal peak inspiratory and expiratory flow and eosinophil cationic proteins (ECPs) were tested for in nasal secretions.	11 participants Age range 18+	FastBlock Allergy Relief had significant effects in reducing symptoms of sneezing and itchy eyes (p<0.01) due to grass pollen allergy. Results indicate that use of FastBlock Allergy Relief can help to alleviate symptoms of hay fever.
Double-blind placebo-controlled trial of cellulose powder as a remedy for persistent allergic rhinitis, by nasal provocation with Der p1 and Der f1 (dust mite allergen) Presented as a Poster EAACI, Gothenburg June 2007 <i>Current Medical Research and Opinion</i> ; Vol 23; No 10; 2007, 2423-2431	Emberlin, Lewis. A trial of FastBlock Allergy Relief by nasal provocation tests with Der p1 and Der f1. Baseline measurements were taken before the challenge of 0.01µg dose containing 5µg of Der p1 and 5µg of Der f1 per g. Measurements taken at 15 minute intervals for the first hour, then at 30 minutes until 4 hours, then at 6 hours and at 24 hours.	15 participants Age range 18+	There was a significant difference in the results for sneezing, itchy nose, runny nose and ECPs in nasal secretions. The peak nasal expiratory and inspiratory flow were also significantly different but there was considerable variation. There was no significant difference between the other symptoms and there were no adverse reactions. FastBlock Allergy Relief can have significant effects in reducing some symptoms of persistent rhinitis due to house dust mite allergy.
Efficacy and safety of medical device FastBlock Allergy Relief in prevention and treatment of persistent allergic rhinitis in adults and children Presented at Moscow XVI Congress for Man and Drugs April 06-10, 2009	Zakharzhevskaya, Sidorenko, Treskunov and Karaulov at the Sechenov Medical Academy, Moscow. The paper describes the findings of an open, non-comparative clinical study of efficacy and safety of FastBlock Allergy Relief in prevention and treatment of persistent allergic rhinitis (AR). Participants were administered FastBlock Allergy Relief 3 times per day over the course of 4-weeks.	48 total participants 25 adults, 23 children Age range 2–62	The severity of AR symptoms and the tolerability of the product were assessed during each visit to the investigator. The results showed that FastBlock Allergy Relief reduces the severity of AR symptoms in the first week of treatment and overall there was significant improvement in symptom reduction by week 4. A twofold improvement in the quality of life of the AR patients was recorded. The researchers concluded that FastBlock Allergy Relief is an effective and safe method of prevention and treatment of allergic rhinitis both in adults and children.
A nasally applied cellulose powder in seasonal allergic rhinitis (SAR) in children and adolescents; reduction of symptoms and relation to pollen load Presented EAACI London 2010 Published <i>Paediatric Allergy and Immunology</i> 22 (2011) 594-599 (2011)	Åberg and Benson. Conducted at the Queen Silvia Children's Hospital, Gothenburg, Sweden in 2009. A double-blind, placebo-controlled trial. All participants were on daily oral antihistamine appropriate for their age and FastBlock Allergy Relief 3 times daily. Reporting of symptoms and reminders were done by SMS.	53 participants Age range 8–18	There was a significant reduction in total symptom scores for the nose (p=0.033) and specifically a runny nose (p=0.0017). There was a tendency for all symptoms scores to be lower for the active group. During low or moderate pollen count there is a significant reduction in total nasal symptom and runny nose along with sneezing severity.
Intranasal Inert Cellulose Powder in Prevention and Management of Seasonal Allergic Rhinitis (SAR) in children Presented as a Poster in EAACI London 2010.	Geppé, Snegotskaya, Kolosova, Konopelko. The study took place at the Clinic of Child Diseases at the I.M Sechenov Moscow Medical Academy in 2009. An open comparative, randomized study. Participants divided into 4 groups depending on their current treatment. Each group received a different medication. Group 1: FastBlock Allergy Relief twice a day. Group 2: Montelukast 5 mg once a day. Group 3: 2 doses of 50 mg Sodium Cromoglicate. Group 4: Budesonide 50 mg 3-4 times a day.	50 participants Age range 8–18	26 patients demonstrated positive results from the very first application of FastBlock Allergy Relief. After 6 weeks treatment, Group 1 demonstrated a significant decrease (all p<0.001) of all SAR symptoms: runny nose, sneezing, nasal blockage, nasal itching, eye itching, nasopharyngeal itching. Children using FastBlock Allergy Relief decreased frequency of use of antihistamine, decongestants and topical steroids. FastBlock Allergy Relief has minimal side effects and is appropriate for children.

STUDY	DESCRIPTION	PARTICIPANTS	RESULTS
Open non-comparative study to evaluate the effectiveness of FastBlock Allergy Relief preparation for patients with allergic rhinitis Published in <i>Russian Allergy Journal</i> in No. 2 (March–April 2011)	Chief clinical physician, N.I. Ilina. The study was conducted at the Russian Federal Medical Biological Agency. An open study over 3 months to determine the effectiveness of FastBlock Allergy Relief at treating allergic rhinitis by nasal provocation test with significantly causative aeroallergens.	30 participants 18 female, 12 male Average age 28.5	The therapy using FastBlock Allergy Relief was found to be effective in 28 (99.6%) of the patients. There was also a significant decrease in nasal reactivity due to a causative allergen. The best results were obtained in patients with isolated dust sensitivity and a mild period of rhinitis. No participants showed any adverse reactions.
Nasal mucociliary clearance and mucoadhesion of hydroxypropylmethylcellulose of powder used for alleviation of allergic rhinitis Presented as a Poster at EAACI, London 2010 Published in <i>Natural Science</i> Vol.2 No.2, 79-84 (2010)	Diethart of School of Human and Health Sciences, Swansea University, Emberlin of National Pollen and Aerobiology Research Unit, University of Worcester and R. Lewis, Worcestershire Royal Hospital. 12 healthy volunteers were tested at the end of the grass pollen season in 2008. The mucociliary clearance time was tested in the absence and then presence of HPMC using a modified Andersen saccharine test.	12 participants 9 female, 3 male Average female age 32.8 Average male age 37.0	When FastBlock Allergy Relief was applied to the nostril, it significantly increased the mucociliary clearance time (MCT). The average MCT at baseline was 11.14 minutes. This significantly increased to 35–45 minutes when 10 mg of FastBlock Allergy Relief was applied to the nostril prior to the test ($p < 0.0005$). Application of 20 mg resulted in an average MCT of 50.37. This increase in MCT was statistically significant when compared to baseline and 10 mg FastBlock Allergy Relief. FastBlock Allergy Relief reduced the mesh spacing of the mucus to form a barrier slowing down MCT.
Efficacy of cellulose powder as part of a complex therapy in patients with intermittent allergic rhinitis <i>Russian Allergy Journal</i> , May 2011	Penechko, Sizyakina. Participants were divided into 2 groups. Group one received standard therapy (second generation cetirizine antihistamine, sorbents and topical glucocorticosteroids). Group two received FastBlock Allergy Relief three times a day in addition to the basic therapy. The observation period was 4 weeks with the patients visiting the clinic once a week.	30 participants Age range 18–33	Group one only saw significant improvement in symptoms such as runny and stuffy nose. Participants in group two showed improvement by the end of the first week and at the end of the fourth week there was a statistically significant reduction in symptoms: runny nose, sneezing, itchy nose and stuffy nose. Analysis of the questionnaire showed significant improvement in quality of life for the second group. Comparative analysis demonstrated that FastBlock Allergy Relief leads to faster alleviation of symptoms and improves quality of life.
A Nasally Applied Cellulose Powder in Seasonal Allergic Rhinitis in Adults with Grass Pollen Allergy: A double-Blind, Randomized, Placebo-Controlled, Parallel-Group study <i>INT Arch Allergy Immunol</i> 2014;163:313 -318	Åberg, Ospanova, Nikitin, Emberlin and Dahl. The study was performed at the University Clinics of Kharkov and Dnepropetrovsk, Ukraine, in May 2013. The patients were randomly assigned FastBlock Allergy Relief or placebo and given identical devices to be puffed in each nostril 3 times daily. Three times a day the patients were reminded by SMS to take their nasal puffs and were asked to confirm the intake with a response SMS. In the evening, they were asked about the severity of symptoms during the day from the nose, eyes, and lower airways and to answer with a figure from 1 (no symptoms) to 6 (strong symptoms).	108 participants Age range 18–40	Significant reductions were found in severity scores for sneezing, runny nose, stuffy nose, and symptoms from eyes and lower airways, both separately and together. Reflective opinion of effect and guess on treatment at follow-up visits confirmed a high efficacy. There was a significant difference between the global opinions of the two groups. 87.1% of the active participants found the product had a good effect. The product provided significant protection against all seasonal allergic rhinitis symptoms. There were no severe adverse events.

Ingredients:

Cellulose powder of vegetable origin and peppermint powder

Uses:

Temporarily relieves mild allergic symptoms of hay fever and allergic rhinitis: sneezing, itchy nasal passages, runny nose, and sinus congestion.

Directions:

Gently blow your nose. While inhaling normally, place tip of dispenser just inside the nostril and firmly squeeze the bottle to distribute powder. Repeat in the other nostril. The smell of peppermint indicates the powder has been delivered.

Adults and Children (18 months and older):

Spray 1–2 times in each nostril, up to 3 times per day.
Children under 18 months old: Consult a physician.

Additional Tips:

Use at least 3 times daily or as often as required. Best results are achieved by maintaining a constant layer of powder across the nasal mucosa. This creates a barrier that prevents allergens from coming into contact with the sensitive membrane within the nasal passage, thereby preventing the symptoms of an allergic reaction.



If you choose to use FastBlock Allergy Relief only when you know you will be encountering allergens (such as mowing the lawn or visiting a home with pets), apply FastBlock Allergy Relief beforehand to establish the barrier.

If you blow your nose, use FastBlock Allergy Relief again to re-establish the barrier.



FASTBLOCK™ ALLERGY RELIEF TESTIMONIALS

Readers of **Woman** Magazine who suffer from hay fever tried FastBlock™ Allergy Relief. Here are results of 190 returned questionnaires.



“Previously, I had used both eye drops and nasal sprays for my hay fever. I had to give those up due to breast feeding and switched to FastBlock Allergy Relief since it is drug-free. I only used it in the morning, and it relieved my symptoms a lot. I would recommend it!”

— J. A.

“A good product which works quickly with no side effects”

“Works instantly, wasn't hard to use”

“I have hay fever during the summer months, and in the past have tried all the usual antihistamines. I tried FastBlock Allergy Relief—it worked within a few minutes! I felt it was better than other allergy solutions because of the speed of relief, how long the relief lasted and it is drug-free! I don't feel like I'm constantly taking drugs anymore. It offered considerably more relief than everything I've taken in the past, and I will be recommending it!” — Wendy P.

“I have suffered with allergies for a long time in the summer, and found it manageable by using a general antihistamine (which didn't really do much). FastBlock Allergy Relief worked quite fast—I felt relief of symptoms within about 15 minutes, and found it better because it worked more quickly than a pill.” — C.O.

“It worked better than any other product I have bought for my hay fever”

“I have suffered from hay fever since the age of ten, along with allergies to animals and dust. Once I started using FastBlock Allergy Relief for a while, the effect was more or less instant and gave me complete relief from symptoms. I use it about every other day, or whenever the need arises. Previously I used prescribed nasal sprays and general antihistamines, but I feel Fastblock Allergy Relief is easier to use, and I like that I don't have to remember if I've had a dose or not, like I needed to recall with the tablets. Now that I'm breastfeeding, the fact that it is drug-free is a definite advantage!”

— Sara R.

“I liked that FastBlock Allergy Relief was easy to use, and I felt quicker relief of symptoms. I recommended it to my father, and he has also used it with great results.” — Sally B.

“I would highly recommend. Great relief all day long from symptoms”

Terry Naturally
EuroPharma.

Our Passion is Your Health!

Manufactured Exclusively for:
EuroPharma, Inc., Green Bay, WI 54311
866-807-2731 | FastBlockAllergy.com